

Careline™ Combi Bag

Non-Sterile

Intended use:

Urine drainage bag with scale for collection and monitoring of urine output. The device is not intended to have a measurement function. The scale is approximate and for indication purposes only. The product is intended for single use only.



Product range:

Careline™ Combi Bag with Non-Return Valve (NRV)

REF*	Capacity, ml	Tube length, cm	Bottom outlet / Cross valve	Cut to fit	Textile on backside
32 450 0-xx	1,500	90	Yes	Yes	Yes

* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains DEHP? No
Tube Polyvinyl chloride
Connector Polyvinyl chloride
Adjustable connector Polypropylene, High density polyethylene
Cap Low density polyethylene, High density polyethylene
Cross valve Polypropylene, High density polyethylene
Film Polyvinyl chloride

Certifications:

CE-marking? Yes
Product classification Class I (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN)58922

Sterilization:

Sterile? No
Shelf life 5 years

Packaging specifications:

REF	Shipping carton				MOQ, pcs
	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Volume m ³	
32 450 0-xx	30	375 x 280 x 81	2.4	0.009	3,840

Storage and handling:

