

Suction Catheters











ProFlo™ and ProFlo™ Plain with Vacutip

Intended use:

Suction catheters are intended for oro-, nasopharyngeal and tracheobronchial suctioning. Must be attached to a jar/canister via tube before being connected to vacuum source. The products are intended for single use only.



Product range:

	Standard Material (Shore A 78) REF*	Description	CH	O.D. mm	Effective Length mm	Color Coding
ProFlo™, Straight tip with Vacutip, Metric cm-marking on catheter						
	2541005101xx	ProFlo™, Straight, Extended length, Metric	10	3.3	570	Black 
	2541006101xx	ProFlo™, Straight, Extended length, Metric	12	4.0	570	White 
	2541007101xx	ProFlo™, Straight, Extended length, Metric	14	4.7	570	Green 
	2541008101xx	ProFlo™, Straight, Extended length, Metric	16	5.3	570	Orange 
ProFlo™ Plain, Straight tip with Vacutip, Metric cm-marking on catheter						
	2541001101xx	ProFlo™ Plain, Straight, Extended length, Metric	10	3.3	570	Black 
	2541002101xx	ProFlo™ Plain, Straight, Extended length, Metric	12	4.0	570	White 
	2541003101xx	ProFlo™ Plain, Straight, Extended length, Metric	14	4.7	570	Green 
	2541004101xx	ProFlo™ Plain, Straight, Extended length, Metric	16	5.3	570	Orange 

* Here noted as xx, reference codes contain two additional digits (10-50) representing a language code. Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
 Contains DEHP? No
 Tube Polyvinyl chloride
 Connector Polyvinyl chloride, Styrene-acrylonitrile copolymer

Certifications:

CE-marking? Yes (CE 0086)
 Product classification Class IIA (EU)
 Quality system certified to ISO 13485
 Global Medical Device Nomenclature (GMDN)10749

Sterilization:

Sterile? Yes
 Sterilized by Ethylene oxide
 Shelf life 5 years

Suction Catheters

ProFlo™ and ProFlo™ Plain with Vacutip



Packaging specifications:

REF	Shelf box		Shipping carton				MOQ, pcs
	Qty, pcs	Dimensions, mm (LxWxD)	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt. kg	Volume m ³	
ProFlo™, Straight tip with Vacutip							
2541005101xx	75	145 x 124 x 734	450	391 x 304 x 743	7.2	0.088	1,350
2541006101xx	75	145 x 124 x 734	450	391 x 304 x 743	7.8	0.088	1,350
2541007101xx	75	145 x 124 x 734	450	391 x 304 x 743	8.5	0.088	1,350
2541008101xx	75	145 x 124 x 734	450	391 x 304 x 743	9.2	0.088	1,350
ProFlo™ Plain, Straight tip with Vacutip							
2541001101xx	75	145 x 124 x 734	450	391 x 304 x 743	7.2	0.088	900
2541002101xx	75	145 x 124 x 734	450	391 x 304 x 743	7.8	0.088	900
2541003101xx	75	145 x 124 x 734	450	391 x 304 x 743	8.5	0.088	900
2541004101xx	75	145 x 124 x 734	450	391 x 304 x 743	9.2	0.088	900

Storage and handling:



ConvaTec Limited

First Avenue
Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK



ConvaTec Limited
First Avenue, Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK
www.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.

All trade marks and logos are the property of a ConvaTec group company.

© 2019 ConvaTec Inc. AP-020199-MM, ver. MAY 2019